

INTRANASAL SPLINTS





PosiSep[™] and PosiSep[™] X are indicated for use as intranasal splints intended to minimize bleeding and edema and to prevent adhesions between the septum and the nasal cavity. They are placed in the nasal cavity after surgery or trauma. The PosiSep™ line is constructed of a patient-comfortable sponge which is manufactured from naturally occurring Chitosan polymers.

- Dissolves naturally or can be easily removed with • gentle irrigation and aspiration
- Biocompatible made of Chitosan-based polymers •
- Minimizes bleeding and edema after surgery •
- Foam easily conforms to nasal anatomy •
- Easy to apply no mixing or additional prep • time necessarv
- No special storage conditions required .

PosiSep[™] and PosiSep[™] X are manufactured from naturally occurring Chitosan polymers.

What is Chitosan?

Chitosan (ky-toe-san) comes from a substance known as chitin, an amino-polysaccharide removed from the powdered shells of crustaceans.



Two Different Platforms for Your Individual Needs

PosiSep X



When PosiSep[™] is hvdrated, its smooth gel action separates sinus mucosal tissue.

When PosiSep[™] X is hvdrated. it expands to fill the space preventing post-op adhesions.

Each PosiSep[™] is available in two sizes. See ordering information below.

ORDERING INFORMATION

Part Number:	Description:
9210574	PosiSep™ Intranasal Splint 0.6" x 2.0" 5pk (US)
9210554	PosiSep [™] Intranasal Splint 0.5" x 1.5" 5pk (US)
9210584	PosiSep [™] X Intranasal Splint 0.6" x 2.0" 5pk (US)
9210564	PosiSep [™] X Intranasal Splint 0.5" x 1.5" 5pk (US)

How Supplied:

PosiSep[™] and PosiSep[™] X are packaged five (5) sterile, single use units per carton. One (1) sterile unit contains one (1) Intranasal Splint. PosiSep[™] and PosiSep[™] X are terminally sterilized by gamma irradiation.

CAUTION: Federal (U.S.) law restricts this device to sale by or on the order of a physician. PosiSep[™] is contraindicated for patients who have known allergies to shellfish. For detailed information regarding indications for use, warnings and precautions, see Instructions for Use.

To place an order or for more information on PosiSep[™] or our other medical technologies, please contact a sales representative by calling 866-612-2568.



5000 Township Parkway St. Paul. MN 55110 Tel 866-612-2568