



# ENT Chairs:

Product Range





|                                      |                   | OP-S4              | OP-S6              | OP-S7                       | OP-S8               | OP-S10             | OP-S14             |
|--------------------------------------|-------------------|--------------------|--------------------|-----------------------------|---------------------|--------------------|--------------------|
| <b>Up/down movement</b>              |                   | Electrical (motor) | Electrical (motor) | Electrical (motor)          | Electrical (motor)  | Electrical (motor) | Electrical (motor) |
| <b>Reclining movement</b>            |                   | Electrical (motor) | Manual (gas)       | Electrical (motor)          | Electrical (motor)  | Electrical (motor) | Electrical (motor) |
| <b>Movement control</b>              | <b>Foot Pedal</b> | Yes                | Yes                | Yes                         | Yes                 | Optional           | Optional           |
|                                      | <b>Keyboard</b>   | --                 | --                 | 2                           | 2                   | 4                  | 4                  |
| <b>Memories</b>                      |                   | --                 | --                 | 6                           | 6                   | 4 / 16             | 4 / 16             |
| <b>Rotation Brake</b>                |                   | Manual             | Manual             | Manual (both sides)         | Manual (both sides) | Electro-magnetic   | Electro-magnetic   |
| <b>Maximum load</b>                  |                   | 150 Kg             | 150 Kg             | 150 Kg                      | 200 Kg              | 200 Kg             | 200 Kg             |
| <b>Fold down armrest</b>             |                   | Yes                | Yes (hide away)    | Yes (hide away)             | Yes (hide away)     | Yes (hide away)    | Yes (hide away)    |
| <b>Rotation</b>                      |                   | 330°               | 360° continuous    | 380° (190° left/190° right) | 300°                | 300°               | 300°               |
| <b>Trendelenburg Position</b>        |                   | Yes                | Yes                | Yes                         | Yes                 | Yes                | Yes                |
| <b>Synchronized sliding backrest</b> |                   | --                 | --                 | --                          | --                  | Yes                | Yes                |
| <b>80° Backrest position</b>         |                   | --                 | --                 | --                          | Yes                 | Yes                | Yes                |
| <b>TFT Monitor</b>                   |                   | --                 | --                 | --                          | --                  | --                 | Optional           |
| <b>Integrated camera</b>             |                   | --                 | --                 | --                          | --                  | --                 | Optional           |
| <b>LED Light source</b>              |                   | --                 | --                 | --                          | --                  | --                 | Optional           |
| <b>Examination Spotlight</b>         |                   | --                 | --                 | --                          | --                  | --                 | Optional           |

# Examination Chairs